

# Test Kit Instructions

October 14, 2020

## **PerkinElmer: AuroFlow AQ Ochratoxin Strip Test** **Using QuickStar Horizon Strip Reader**

### **Table of Contents**

<b>GENERAL INFORMATION.....</b>	<b>2</b>
<b>PREPARATION OF TESTING MATERIALS AND EQUIPMENT .....</b>	<b>3</b>
<b>SAMPLE PREPARATION AND EXTRACTION PROCEDURES.....</b>	<b>4</b>
<b>TESTING PROCEDURE.....</b>	<b>4</b>
<b>REPORTING AND CERTIFYING TEST RESULTS.....</b>	<b>7</b>
<b>STORAGE CONDITIONS AND PRECAUTIONS.....</b>	<b>7</b>
<b>EQUIPMENT AND SUPPLIES .....</b>	<b>8</b>
<b>REVISION HISTORY .....</b>	<b>9</b>

## GENERAL INFORMATION

The AuroFlow AQ Ochratoxin Strip Test is a quantitative rapid lateral flow assay (strip test) kit designed to quantitatively determine the amount of Ochratoxin (OTA) residues in wheat. This kit utilizes an environmentally-friendly aqueous extraction method. Ochratoxin present in the sample extraction interacts with the strip (device) components and influences color intensity on the test and control lines. Developed test strips are inserted into an electronic strip reader, which reports Ochratoxin contamination of the sample in parts per billion (ppb). This test is designed for rapid field or reference laboratory settings.

Please read all instructions thoroughly prior to performing the test.

The instructions presented in this document cover only the procedure for performing the analytical test for official inspections. For questions regarding this procedure, contact Dr. Ajit Ghosh of the Technology and Science Division by phone at 816-891-0417 or email at [Ajit.K.Ghosh@usda.gov](mailto:Ajit.K.Ghosh@usda.gov).

Refer to the Mycotoxin Handbook for information on the use of this test kit for official inspections including sampling, general sample preparation, reporting and certification of test results, laboratory safety, and hazardous waste management. For questions regarding these policies and/or instructions, contact Patrick McCluskey (816-659-8403 or [Patrick.J.McCluskey@usda.gov](mailto:Patrick.J.McCluskey@usda.gov)).

### Approved Test Kit Information

Kit Test Vendor:	PerkinElmer Inc. (1-512-707-8993)
Test Kit Name:	AuroFlow AQ Ochratoxin A Strip Test
Product Number:	FOOD-1417-01
Effective Date of Instructions:	10/14/2020
Conformance Range:	5.0 – 100 ppb
Number of Analyses to Cover Conformance Range:	1
Type of Service:	Quantitative
Approved Commodities:	Wheat (including whole grain wheat flour only)
Extraction Method:	Shake 50 grams of sample with 200 mL AQ Extraction Buffer for 1 minute and 30 seconds.
Test Format:	Lateral Flow Strip
Detection Method:	QuickStar Horizon Strip Reader (Catalog # FOOD-6006-01)

## PREPARATION OF TESTING MATERIALS AND EQUIPMENT


### 1. Equilibration of Kit Contents

- a. Remove the kit and all its accompanying contents from 4 °C storage. Allow the kit to rest for 1 hour at room temperature (18 – 30 °C) before performing any determinations.
- b. Do not use the kit outside an 18 – 30 °C (65 – 86 °F) temperature range.

### 2. Preparation of AQ Extraction Buffer:

- a. AQ Extraction Buffer must be prepared prior to sample extraction.
- b. To prepare the AQ Extraction Buffer, add 550 mL of distilled or deionized water using a 1 L graduated cylinder and 300 µL of AQ Add-In Solution using a 1000 µL pipettor into a clean 1 L bottle.
- c. Transfer the entire contents of one (1) Large AQ Powder Pack (Catalog # FOOD-4005-200) into the bottle with a wide-mouth funnel and tap to dislodge solids. Manually shake for 1 minute to dissolve completely. Label it with the date and the analyst's name. Store at room temperature for up to 3 months.
  1. Alternatively: Add a magnetic stir bar and stir at high speed with a magnetic stir plate.

### 3. Preparation of Equipment

- a. Turn on orbital shaker and confirm calibration status.
- b. Plug the reader and printer into their respective power supplies. Power on the reader. Connect the printer to the reader via a USB cord.
- c. Calibrate the strip reader before performing any extractions:
  1. QuickStar Horizon Strip Reader
    - a. Power on the reader.
    - b. Select user or enter user information.
    - c. Identify the QR Calibration Code matching the kit lot #.
    - d. From the main menu, navigate to settings via the gear icon .
    - e. Select "Lot ID".
    - f. Scan the supplied QR code. The lot-specific calibration information will automatically be saved to the matching Lot ID #.

## SAMPLE PREPARATION AND EXTRACTION PROCEDURES

1. Extraction procedure for **wheat** and **whole grain wheat flour**. The sample to be tested should be collected and prepared according to accepted sampling techniques (see **Mycotoxin Handbook**).
  - a. Transfer  $50 \pm 0.2$  grams of ground sample into a Large Extraction Bag (Catalog # FOOD-4009-200)
  - b. Transfer 200 mL of AQ Extraction Buffer into container using a 250 mL graduated cylinder.
  - c. Seal the Large Extraction Bag and briefly shake to moisten the sample and break the lump.
    1. Shake the sample at 250 rpm with an orbital shaker for 1 minute and 30 seconds.
    2. To clarify the extract, prepare the filtration unit: Fold one filter (Catalog # FOOD-4010-01) and place it inside of a filter funnel. Then place the funnel into a clean conical tube.
    3. Using a serological or transfer pipet, transfer approximately 10 mL of liquid extract into the assembled filter funnel. Avoid transferring solids or foam.
    4. Allow the sample to filter by gravity until approximately 2 – 3 mL of filtrate has been collected (approximately 7 minutes). This filtrate is the “**Sample Extract**”.
    5. Remove the funnel and filter from the conical tube and set aside. Residual solids and waste on the filter can be discarded.
    6. Cap the conical tube, and then briefly vortex the contents to mix.
    7. The “**Sample Extract**” must be tested within 30 minutes, as described below.
    8. Follow local and organizational guidance for disposing of Ochratoxin - contaminated samples.

## TESTING PROCEDURE

1. Reader Setup
  - a. Select the appropriate program on the QuickStar Horizon Strip Reader and select the correct lot ID #. Enter any desired sample information.
2. Preparation of **Sample Mix**


- a. Quantitation Range (5-100 ppb): All samples must be tested using the following procedure:
  1. Using the 20 - 200  $\mu\text{L}$  variable pipettor, transfer 150  $\mu\text{L}$  of Running Buffer to a capped tube.
  2. Using the 20 - 200  $\mu\text{L}$  variable pipettor, transfer 150  $\mu\text{L}$  of the **Sample Extract** into the same tube containing Running Buffer. Pipet up and down 3 times, and then vortex on high for 5 seconds to mix. This is the **"Sample Mix"**.
  3. Test the **Sample Mix** within 30 minutes, as described in "Sample Analysis".

### 3. Sample Analysis

- a. Using the 20 - 200  $\mu\text{L}$  variable pipettor, transfer 150  $\mu\text{L}$  of the **Sample Mix** to an appropriate Capless Tube in a tube rack at room temperature.
- b. Add a new test strip to the Capless Tube with the arrows pointing down and allow the strip to develop for 5 minutes.
- c. Remove the strip from the Capless Tube immediately at 5 minutes. Gently touch the end of the test strip onto an absorbent paper towel or Kimwipe.
- d. Visually inspect strips immediately (**Figure 1**), as described below. Valid strips must be analyzed within 1 minute of removal. Do not read invalid strips.

### 4. Visual Inspection





- a. The test strip is INVALID if any of the following are observed:
  1. Test (T) or Control (C) line is uneven or "broken".
  2. Control (C) line is absent.
  3. Gold (pink conjugate) does not flow past T and C lines.



*Figure 1. Test Strip Orientation and Test Interpretation*

- b. Do not read INVALID test strips in reader.
- c. If test strip is invalid, prepare fresh **Sample Mix** and re-test with a new strip.
- d. If test strip is valid (no observed defects and even line development), proceed to “Interpretation” below to read strip.

## 5. Interpretation

- a. Only use readers calibrated with a lot-specific QR code matching the strip lot.
- b. QuickStar Horizon Strip Reader.
  - 1. Insert a valid strip into the cassette, with arrows inserted first, and with “#1” facing up. Select “Run Test”. Select correct Lot ID #.
  - 2. If desired, enter sample ID.
  - 3. The reader should automatically record a result.
    - a. If the reader does not automatically read the strip, select the “” icon to read the strip manually.
  - 4. Verify that the reader displays both control and test lines in the window, otherwise the test is invalid. The strip may be re-read if necessary using the “” icon.
    - a. Invalid results require repeating the analysis from the “Sample Mix Preparation” step.
  - 5. Print the result by pressing the printer icon “”.
- c. If the QuickStar Horizon Strip Reader reports an “Invalid” result, verify that the test strip is properly inserted into the cassette and reader. If the reader again reports an “Invalid” result, perform a calibration verification as follows:
  - 1. In the top left Settings “” menu select “Diagnostic Check”.
  - 2. Following the on-screen instructions insert the white crosshair cartridge into the cassette slot, rounded end down. Make sure the cartridge is firmly seated all the way at the bottom of the slot.
  - 3. Press “Start” and allow the reader to complete the diagnostic.
  - 4. Confirm that the date of the last diagnostic has changed to the current date and the status is displayed as “Valid”.
  - 5. If the reader continues to report “Invalid” results, contact [Bioo.Support@PerkinElmer.com](mailto:Bioo.Support@PerkinElmer.com) for support.

## REPORTING AND CERTIFYING TEST RESULTS

Refer to the Mycotoxin Handbook for reporting and certification of test results. For questions regarding these instructions, contact Patrick McCluskey (816-659-8403 or [Patrick.J.McCluskey@usda.gov](mailto:Patrick.J.McCluskey@usda.gov)).

## STORAGE CONDITIONS AND PRECAUTIONS

### 1. Storage Conditions:

- a. Store the test strips refrigerated at 4 °C in the closed original container.
- b. Store the Running Buffer bottles refrigerated at 4 °C.
- c. Store Large AQ Powder Packs and AQ Add-in Solution with the kit at 4 °C until preparation of AQ Extraction Buffer.
- d. Store prepared AQ Extraction Buffer at room temperature (18 – 30 °C) for up to 3 months. Discard after 3 months or if the solution becomes turbid before 3 months.

### 2. Precautions:

- a. Do not leave test strip canister open for prolonged periods of time; remove strips as needed and re-close the canister.
- b. Do not perform determinations without equilibrating the test kit to room temperature.
- c. Do not perform more than 4 determinations simultaneously.
- d. Only use Running Buffer with test strips from the matching kit lot.
- e. AQ Extraction Buffer (not lot specific) is non-hazardous and may be disposed as normal waste.
- f. Do not open the Large AQ Powder Packs until ready to use.
- g. Do not use test kits or components beyond the expiration date.
- h. The QuickStar Horizon Strip Reader cassette may accumulate residue after extended use, in which case it may be wiped clean with distilled water applied to a Kimwipe. This residue build-up may be reduced by gently touching the freshly-developed strip on a paper towel or Kimwipe before inserting into the cassette.

## EQUIPMENT AND SUPPLIES

1. Required materials provided in the AuroFlow AQ Ochratoxin Strip Test (Catalog # FOOD-1417-01) kit:
  - a. 25 Test Strips
  - b. 1 bottle of Running Buffer
  - c. 25 Filters
  - d. 25 Capless Tubes
  - e. 25 Capped Tubes
  - f. 25 Transfer Pipets, 3 mL
  - g. 1 Calibration Code (lot-specific)
2. Required materials provided in the AuroFlow AQ Supplemental Extraction Reagents (Catalog # FOOD-141301-200) kit:
  - a. 25 Large Extraction Bags
  - b. 9 Large AQ Powder Packs
  - c. 3 tubes of AQ Add-in Solution
3. Required materials, but not provided:
  - a. Strip Reader
    1. QuickStar Horizon Strip Reader (Catalog # FOOD-6006-01)
  - b. Recommended orbital shaker setup:
    1. Orbital shaker (VWR # 0027-142)
    2. Shaker platform (VWR # 97003-590)
    3. Stainless Steel Clamps (VWR # 14215-238)
    4. 1000 mL Plastic Beaker (VWR # 76266-414)
  - c. Microtube rack
  - d. Distilled or deionized water
  - e. Timer or watch
  - f. Single-channel pipettor(s)



1. 20 - 200  $\mu$ L variable pipettor: used for 150  $\mu$ L and 300  $\mu$ L (two times 150  $\mu$ L for AQ Add-In Solution) volume transfers (VWR # 89079-970)
2. 200  $\mu$ L pipet tips (VWR # 89136-161)
- g. Pipet controller
- h. 10 mL graduated serological pipets (VWR # 89130-898)
- i. Scale capable of measuring 50 grams – see Mycotoxin Handbook
- j. Grain mill / grinder – see Mycotoxin Handbook
- k. 250 mL graduated cylinder(VWR #65000-008)
- l. 1 L graduated cylinder (VWR #65000-012)
- m. Storage bottle(s) for AQ Extraction Buffer (glass or plastic, VWR # 16157-282)
- n. Vortex mixer (VWR # 10153-838)
- o. Filter funnel (VWR # 414004-288)
- p. Conical tubes, 50 mL (VWR # 89174-474)
- q. Wide-mouth funnel (VWR # 414004-272)
- r. AuroFlow AQ Supplemental Extraction Reagents (Catalog # FOOD-141301-200)
4. Optional materials not provided:
  - a. QuickStar Horizon Printer (Catalog # FOOD-6007-01)

## REVISION HISTORY

Effective: 10/14/2020